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COMPEX Wireless: your personal coach online

COMPEX WIRELESS

## CREATE YOUR ACCOUNT

test.complex.sp8@gmail.com

Stepan G

Ukraine

•••••

•••••

male  female

1984

Yes, please keep me updated on Complex news

I accept terms and conditions\*

**REGISTER**

\*indicates required field

COMPEX Wireless your personal coach online

WELCOME

My Objectives

- Improve quadriceps strength (Explosive strength)
- Increase your muscle volume (Muscle building)
- My planning
- Optimize my overcompensation (Posture)

My Programmes

Active history

- Explosive strength 1
- Part 2
- Part 3
- Muscle building 1
- Posture

|                      | Monday | Tuesday | Wednesday | Thursday | Friday | Saturday | Sunday |
|----------------------|--------|---------|-----------|----------|--------|----------|--------|
| Explosive strength 1 |        |         |           |          |        |          |        |
| Part 2               |        |         |           |          |        |          |        |
| Part 3               |        |         |           |          |        |          |        |
| Muscle building 1    |        |         |           |          |        |          |        |
| Posture              |        |         |           |          |        |          |        |

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The image shows a web browser window displaying the registration confirmation page for the CompeX website. The browser's address bar shows the URL: `https://www.compexwireless.com/en_EU/auth/confirmation-sent`. The page features a dark blue header with the CompeX logo and navigation links: HOME, PRODUCT INFO, CHAMPIONS STORIES, EMS, FAQ, MEDIA, WIRELESS COACH, and STORE. A banner image shows four CompeX cycling computers (SP 6.0, SP 8.0, Fil 5.0, and another SP 6.0) and a muscular cyclist. A blue circular badge on the right says "European consumer choice". On the left, there are orange navigation buttons: COMPEX STORE, COMPEX PRODUCTS, COMPEX WIRELESS, and WIRELESS COACH. The main content area is dark blue with the text: "Thanks for registering. Confirmation has been sent to your e-mail address."

your personal | X +

os://www.compexwireless.com/en\_EU/auth/confirmation-sent

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COMPEX STORE CHOOSE LANGUAGE SHARE f | t

COMPEX  
gent training partner

DISCOVER THE  
COMPEX WORLD

COMPEX  
STORE

COMPEX  
PRODUCTS

COMPEX  
WIRELESS

WIRELESS  
COACH

HOME PRODUCT INFO CHAMPIONS STORIES EMS FAQ MEDIA WIRELESS COACH STORE

European  
consumer  
choice

Thanks for registering. Confirmation has been sent to your e-mail address.

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The screenshot shows a Gmail interface with a confirmation email from 'Compex' (compexwsnoreply@contentforces.pl) received 13:42 (5 minutes ago). The email content is as follows:

Hello,

This is a confirmation e-mail.

Thank you for registering. To confirm your e-mail login, follow the instructions below.

CONFIRM E-MAIL AND LOGIN

1. Go to the following link:  
[https://www.compexwireless.com/en\\_EU/auth/confirmation/token/ebff9df7-473f-42de-8cef-16e2cc358ed4](https://www.compexwireless.com/en_EU/auth/confirmation/token/ebff9df7-473f-42de-8cef-16e2cc358ed4)
2. Follow the instructions on the web page that opens.

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The image shows a browser window with the URL [https://www.complexwireless.com/en\\_EU/auth/login](https://www.complexwireless.com/en_EU/auth/login). The website header includes the Complex logo, navigation links (HOME, PRODUCT INFO, CHAMPIONS STORIES, EMS, FAQ, MEDIA, WIRELESS COACH, STORE), and utility links (COMPLEX STORE, CHOOSE LANGUAGE, SHARE). A banner features four fitness trackers (SP 6.0, SP 8.0, Fit 5.0, and WIRELESS) and a muscular man. A 'European consumer choice' badge is also present.

On the left sidebar, there are four orange buttons: COMPLEX STORE, COMPLEX PRODUCTS, COMPLEX WIRELESS, and WIRELESS COACH.

The main content area displays the login form with the text "ALREADY REGISTERED? PLEASE LOGIN:". The form includes a text input for the email address (test.compex.sp8@gmail.com), a password input field, a "Remember me" checkbox, a "Password forgotten" link, and a "Sign Up" button. A blue "GO" button is highlighted with a red border.

To the right, a preview of the user dashboard is shown. It features a "WELCOME" message, "My Objectives" (e.g., Improve endurance strength, Increase your muscle volume), "My Programmes" (e.g., Active recovery, Endurance strength 1), and a calendar view for December 2014.

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our personal X +

://www.compexwireless.com/en\_EU/my-settings/chose-device

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COMPEX STORE CHOOSE LANGUAGE SHARE f | t

mpex  
ent training partner

DISCOVER THE  
COMPEX WORLD

COMPEX STORE

COMPEX PRODUCTS

COMPEX WIRELESS

WIRELESS COACH

HOME PRODUCT INFO CHAMPIONS STORIES EMS FAQ MEDIA WIRELESS COACH STORE

European consumer choice

Please select your Compex Wireless device

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Wireless: your personal X +

https://www.compexwireless.com/en\_EU/support/support

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COMPEx STORE CHOOSE LANGUAGE SHARE f | t

HOME PRODUCT INFO CHAMPIONS STORIES EMS FAQ MEDIA WIRELESS COACH STORE

Compex®  
Your intelligent training partner

DISCOVER THE COMPEX WORLD

COMPEX STORE

COMPEX PRODUCTS

COMPEX WIRELESS

WIRELESS COACH

MY DASHBOARD MY SETTINGS QUICK START GUIDE **SUPPORT** LOGOUT

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

The screenshot shows a web browser window with the URL [www.compexwireless.com/en\\_EU/my-settings](http://www.compexwireless.com/en_EU/my-settings). The page features a navigation menu with options: MY DASHBOARD, MY SETTINGS (highlighted with a red box), QUICK START GUIDE, SUPPORT, and LOGOUT. Below the navigation, there are links for MY SETTINGS and RESET PASSWORD. The main content area is titled "WELCOME Stepan G" and contains two sections: "Device settings" and "User's settings".

**Device settings**

|                      |                     |                        |
|----------------------|---------------------|------------------------|
| Device Type          | SP8.0               | <a href="#">change</a> |
| Paired Device S/N    | YRP020157           |                        |
| Language             | English             |                        |
| Last synchronization | 2017-05-13 12:17:46 |                        |

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your personal X +

https://www.complexwireless.com/en\_EU/dashboard#calendar\_

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DISCOVER THE COMPLEX WORLD

COMPLEX STORE

COMPLEX PRODUCTS

COMPLEX WIRELESS

WIRELESS COACH

MY DASHBOARD

MY SETTINGS

QUICK START GUIDE

SUPPORT

LOGOUT

SCHEDULE

CALENDAR

ADD OBJECTIVE

WELCOME

Stephan D

My Objectives

My planning

11-05-2017

I want to add:

Objective

Programme

Task

Note

Ankle Twist Prevention

Save

Friday

Saturday

Sunday

5 +

6 +

7

12 +

13 +

14

+

15 +

16 +

17 +

18 +

19 +

20 +

21

+

22 +

23 +

24 +

25 +

26 +

27 +

28

+

29 +

30 +

31 +

+

+

+

+

# WELCOME

Stepan G

## My Objectives

My planning

Please use the selector to find your objective

Select a sport

Running half marathon

Select a category

Conditioning

Select a body zone

- + Improve my core stabilisation
- + Optimise my overcompensation
- + Improve quadriceps endurance
- + Prepare my half-marathon (3 trainings/week)
- + Prepare my half-marathon (5 trainings/week)

Prepare my half-marathon (3 trainings/week)

Add to my objectives

Objectives

Prepare for a half marathon by combining electrostimulation and running.

Benefits of

Electrostim

the muscle

recovery p

recover fa

Explanation

This objec

The transp

factors wh

Voluntary

these two

suitable fo

their capacity to use oxygen in larger quantities, which contributes to improving endurance performance. Adding a capillarisation session after your long weekly run also enables you to develop the network of microscopic blood vessels (capillaries) found in your muscles to further increase your level of endurance. Post-training recovery will be significantly accelerated by doing a Compex recover session.

The last week before the race the overcompensation effect triggered and improved with several electrostimulation sessions helps you to be ready for the race.

### Add to my objectives

Name\*

Prepare my half-marathon (3

Date\*

11.05.2017

- Start date
- End date

Cancel

Create an Objective

MY DASHBOARD

MY SETTINGS

QUICK START GUIDE

SUPPORT

LOGOUT

SCHEDULE

CALENDAR

ADD OBJECTIVE

Objective has been added successfully

# WELCOME

Stepan G

## My Objectives

- My planning
- Prepare my half-marathon (3 trainings/week)  
End date 14.07.17

| ◀ May 2017 ▶                       |                 |               |          |                                  |          |                 |
|------------------------------------|-----------------|---------------|----------|----------------------------------|----------|-----------------|
| Monday                             | Tuesday         | Wednesday     | Thursday | Friday                           | Saturday | Sunday          |
| 1                                  | 2               | 3             | 4        | 5                                | 6        | 7               |
| 8                                  | 9               | 10            | 11       | 12                               | 13       | 14<br>Endurance |
| 15<br>Running<br>Training recovery | 16<br>Endurance | 17<br>Running | 18       | 19<br>Running<br>Capillarization | 20       | 21<br>Endurance |
| 22<br>Running<br>Training recovery | 23<br>Endurance | 24<br>Running | 25       | 26<br>Running<br>Capillarization | 27       | 28<br>Endurance |
| 29<br>Running<br>Training recovery | 30<br>Endurance | 31<br>Running | 1        | 2<br>Running<br>Capillarization  | 3        | 4<br>Endurance  |
| 5<br>Running<br>Training recovery  | 6<br>Endurance  | 7<br>Running  | 8        | 9<br>Running<br>Capillarization  | 10       | 11<br>Endurance |

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The screenshot shows a web browser window with the URL [www.complexwireless.com/en\\_EU/dashboard](http://www.complexwireless.com/en_EU/dashboard). The page features a navigation menu with options like 'MY DASHBOARD', 'MY SETTINGS', 'QUICK START GUIDE', 'SUPPORT', and 'LOGOUT'. Below this, there are links for 'SCHEDULE', 'CALENDAR', and 'ADD OBJECTIVE'. A message 'Item has been deleted' is visible. The main content area is titled 'WELCOME' and includes a 'My Objectives' section with a list of tasks. A calendar grid is displayed, showing dates from Friday to Sunday. A modal window is open over the calendar, titled '17-05-2017', and contains the text 'I want to add:'. Underneath, there are four options: 'Objective', 'Programme', 'Task', and 'Note'. The 'Objective' option is selected, and a dropdown menu shows 'Get a 6 pack'. A 'Save' button is also visible.

our persona X +

www.complexwireless.com/en\_EU/dashboard

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European consumer choice

MY DASHBOARD MY SETTINGS QUICK START GUIDE SUPPORT LOGOUT

SCHEDULE CALENDAR ADD OBJECTIVE

Item has been deleted

WELCOME

My Objectives

- My planning
- Prepare my half-marathon (trainings/week)  
End date: 14.07.17

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www.compexwireless.com/en\_EU/dashboard

SCHEDULE CALENDAR ADD OBJECTIVE

Program has been added successfully

WELCOME  
Stepan G

My Objectives

- My planning
- Prepare my half-marathon (trainings/week)  
End date 14.07.17

17-05-2017

I want to add:

- Objective
- Programme Ankle Twist Prevention Save
- Task**
- Note**

Monday Tuesday Wednesday Thursday Friday Saturday Sunday

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Open [www.compexwireless.com](http://www.compexwireless.com)

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